



The 3rd Largest Stapling  
Manufacturer in the World

# Ordering Catalogue



Official UK Distributor





**Panther**®  
**Healthcare**

**Established  
in 2002**



**Present  
in over 80  
countries**



**Hubs in North America,  
South America, Europe,  
Middle East,  
and Asia**



**Green Path  
Project**



**Global surgical  
training programme**



**Panther  
Surgical  
Academy**

## Our focuses



### **Innovation**

Developing new  
technology in-line with  
the latest surgical  
stapling demands



### **Quality**

Ensuring the highest  
standards of quality for  
better patient safety




### **Support**

Offering high levels of  
training and expertise  
for our customers



5



8



9



12

For all orders and enquiries  
contact LaproSurge at  
[info@laprosurge.com](mailto:info@laprosurge.com)

## Contents

### Endo Linear Stapling

Endo Linear Cutter Staplers 4  
Endo Linear Cutter Reloads 5

4

4

5

### Open Stapling

Linear Cutter Staplers 8  
Linear Staplers 9  
Circular Staplers 12

8

8

9

12

### Haemorrhoidal Stapling

14

### Hernia Fixation

15

See our product range



# Endo Linear Stapling



## Endo Linear Cutter Staplers

### SMART Plus Fully Powered Handle

| Code | Description                                          | Shaft Length (mm) | Box Quantity |
|------|------------------------------------------------------|-------------------|--------------|
| PEPS | SMART Plus Fully Powered Stapler (Short Handle)      | 65                | 1            |
| PEPM | SMART Plus Fully Powered Stapler (Standard Handle)   | 155               | 1            |
| PEPL | SMART Plus Fully Powered Stapler (Extra Long Handle) | 255               | 1            |

### SMART Plus Half Powered Handle

| Code | Description                                         | Shaft Length (mm) | Box Quantity |
|------|-----------------------------------------------------|-------------------|--------------|
| PECS | SMART Plus Half Powered Stapler (Short Handle)      | 65                | 1            |
| PECM | SMART Plus Half Powered Stapler (Standard Handle)   | 155               | 1            |
| PECL | SMART Plus Half Powered Stapler (Extra Long Handle) | 255               | 1            |

### Endo Manual Linear Cutter Handle

| Code    | Description                                   | Shaft Length (mm) | Box Quantity |
|---------|-----------------------------------------------|-------------------|--------------|
| CEAC-30 | Manual Endo Linear Cutter (Standard Handle)   | 160               | 1            |
| CEAC-45 | Manual Endo Linear Cutter (Short Handle)      | 60                | 1            |
| CEAC-60 | Manual Endo Linear Cutter (Extra Long Handle) | 260               | 1            |

# Endo Linear Stapling

## Endo Linear Cutter Reloads

### P2G Articulating Cartridge with Double Grip Technology (Straight Tip)

| Code       | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Trocar Diameter (mm) | Box Quantity |
|------------|--------------------|-------------------------|---------------------------|----------------------|--------------|
| ● CADF-30T | 30                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADF-30N | 30                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADF-45T | 45                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADF-45N | 45                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADF-45S | 45                 | 4.2                     | 1.8                       | 12                   | 6            |
| ● CADF-45D | 45                 | 4.8                     | 2.0                       | 15                   | 6            |
| ● CADF-45R | 45                 | 5.0                     | 2.2                       | 15                   | 6            |
| ● CADF-60T | 60                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADF-60N | 60                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADF-60S | 60                 | 4.2                     | 1.8                       | 12                   | 6            |
| ● CADF-60D | 60                 | 4.8                     | 2.0                       | 15                   | 6            |
| ● CADF-60R | 60                 | 5.0                     | 2.2                       | 15                   | 6            |




# Endo Linear Stapling



## P2G Articulating Cartridge with Double Grip Technology (Curved Tip)


| Code       | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Trocar Diameter (mm) | Box Quantity |
|------------|--------------------|-------------------------|---------------------------|----------------------|--------------|
| ● CADJ-30T | 30                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADJ-30N | 30                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADJ-45T | 45                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADJ-45N | 45                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADJ-45S | 45                 | 4.2                     | 1.8                       | 12                   | 6            |
| ● CADJ-45D | 45                 | 4.8                     | 2.0                       | 12                   | 6            |
| ● CADJ-45R | 45                 | 5.0                     | 2.2                       | 15                   | 6            |
| ● CADJ-60T | 60                 | 2.5                     | 1.0                       | 12                   | 6            |
| ● CADJ-60N | 60                 | 3.5                     | 1.5                       | 12                   | 6            |
| ● CADJ-60S | 60                 | 4.2                     | 1.8                       | 12                   | 6            |
| ● CADJ-60D | 60                 | 4.8                     | 2.0                       | 12                   | 6            |
| ● CADJ-60R | 60                 | 5.0                     | 2.2                       | 15                   | 6            |



**CADJ-60N**

## 8mm Articulating Reload

| Code       | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Trocar Diameter (mm) | Box Quantity |
|------------|--------------------|-------------------------|---------------------------|----------------------|--------------|
| ● CADL-45T | 45                 | 2.5                     | 1.0                       | 10                   | 6            |



**CADL-45T**

# Endo Linear Stapling

**CADF-45S**



**PEPS**



**PECM**



**CEAC-30**




# Open Stapling



## Linear Cutter Stapler


| Code        | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Box Quantity |
|-------------|--------------------|-------------------------|---------------------------|--------------|
| ● SSAB-60N  | 60                 | 3.8                     | 1.5                       | 1            |
| ● SSAB-60D  | 60                 | 4.8                     | 2.0                       | 1            |
| ● SSAB-80N  | 80                 | 3.8                     | 1.5                       | 1            |
| ● SSAB-80D  | 80                 | 4.8                     | 2.0                       | 1            |
| ● SSAB-100N | 100                | 3.8                     | 1.5                       | 1            |
| ● SSAB-100D | 100                | 4.8                     | 2.0                       | 1            |



**SSAB-60D**

## Linear Cutter Stapler Reloads

| Code        | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Box Quantity |
|-------------|--------------------|-------------------------|---------------------------|--------------|
| ● SADB-60N  | 60                 | 3.8                     | 1.5                       | 3            |
| ● SADB-60D  | 60                 | 4.8                     | 2.0                       | 3            |
| ● SADB-80N  | 80                 | 3.8                     | 1.5                       | 3            |
| ● SADB-80D  | 80                 | 4.8                     | 2.0                       | 3            |
| ● SADB-100N | 100                | 3.8                     | 1.5                       | 3            |
| ● SADB-100D | 100                | 4.8                     | 2.0                       | 3            |




**SADB-60N**

# Open Stapling



## Linear Stapler

| Code            | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Box Quantity |
|-----------------|--------------------|-------------------------|---------------------------|--------------|
| ● FLSLF30 3.5mm | 30                 | 3.5                     | 1.5                       | 1            |
| ● FLSLF30 4.8mm | 30                 | 4.8                     | 2.0                       | 1            |
| ● FLSLF45 3.5mm | 45                 | 3.5                     | 1.5                       | 1            |
| ● FLSLF45 4.8mm | 45                 | 4.8                     | 2.0                       | 1            |
| ● FLSLF60 3.5mm | 60                 | 3.5                     | 1.5                       | 1            |
| ● FLSLF60 4.8mm | 60                 | 4.8                     | 2.0                       | 1            |
| ● FLSLF90 3.5mm | 90                 | 3.5                     | 1.5                       | 1            |
| ● FLSLF90 4.8mm | 90                 | 4.8                     | 2.0                       | 1            |



**FLSLF30 3.5mm**

# Open Stapling



## Linear Stapler Reloads

| Code          | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Box Quantity |
|---------------|--------------------|-------------------------|---------------------------|--------------|
| ● (FLSL)H3.5A | 30                 | 3.5                     | 1.5                       | 3            |
| ● (FLSL)H4.8A | 30                 | 4.8                     | 2.0                       | 3            |
| ● (FLSL)H3.5V | 45                 | 3.5                     | 1.5                       | 3            |
| ● (FLSL)H4.8B | 45                 | 4.8                     | 2.0                       | 3            |
| ● (FLSL)H3.5C | 60                 | 3.5                     | 1.5                       | 3            |
| ● (FLSL)H4.8C | 60                 | 4.8                     | 2.0                       | 3            |
| ● (FLSL)H3.5E | 90                 | 3.5                     | 1.5                       | 3            |
| ● (FLSL)H4.8E | 90                 | 4.8                     | 2.0                       | 3            |

**(FLSL)H3.5A**



# Open Stapling



## Linear Stapler with Knife

| Code       | Description               | Staple Length (mm) | Staple Height Open (mm) | Staple Height Closed (mm) | Box Quantity |
|------------|---------------------------|--------------------|-------------------------|---------------------------|--------------|
| ● CAAB-40  | Linear Stapler with Knife | 40                 | 4.8                     | 2.0                       | 1            |
| ● CAZD-40D | Reloading Unit            | 40                 | 4.8                     | 2.0                       | 1            |



**CAZD-40D**



**CAAB-40**

# Circular Staplers



## Circular Stapler with Fixed Anvil

| Code      | Shaft Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|-----------|-------------------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSLWBE21 | 216               | 21                  | 13                  | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSLWBE24 | 216               | 24                  | 15                  | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSLWBE26 | 216               | 26                  | 17                  | 20                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSLWBE29 | 216               | 29                  | 20                  | 24                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSLWBE32 | 216               | 32                  | 22                  | 30                | 5.0                     | 2.2                       | 0.28             | 1            |
| FCSLWBE34 | 216               | 34                  | 24                  | 32                | 5.0                     | 2.2                       | 0.28             | 1            |

## Circular Stapler with Folding Anvil

| Code    | Shaft Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|---------|-------------------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSME21 | 213               | 21                  | 13                  | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSME24 | 213               | 24                  | 15                  | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSME26 | 213               | 26                  | 17                  | 20                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSME29 | 213               | 29                  | 20                  | 24                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSME32 | 213               | 32                  | 22                  | 30                | 5.0                     | 2.2                       | 0.28             | 1            |

# Circular Staplers



## Circular Stapler with Folding Anvil & Feeding Tube Solution

| Code      | Shaft Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|-----------|-------------------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSSWAF24 | 212               | 24                  | 15                  | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSSWAF26 | 212               | 26                  | 17                  | 20                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSSWAF29 | 212               | 29                  | 20                  | 24                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSSWAF32 | 212               | 32                  | 22                  | 30                | 5.0                     | 2.2                       | 0.28             | 1            |

## XL Circular Stapler with Folding Anvil & Feeding Tube Solution

| Code      | Shaft Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|-----------|-------------------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSLWAF21 | 324               | 21                  | 12.9                | 18                | 4.5                     | 1.8                       | 0.28             | 1            |
| FCSLWAF25 | 324               | 25                  | 17                  | 20                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSLWAF29 | 324               | 29                  | 20                  | 24                | 4.8                     | 2.0                       | 0.28             | 1            |
| FCSLWAF32 | 324               | 32                  | 22                  | 30                | 5.0                     | 2.2                       | 0.28             | 1            |
| FCSLWAF33 | 324               | 33                  | 24                  | 32                | 5.0                     | 2.2                       | 0.28             | 1            |

# Haemorrhoidal Staplers



## Haemorrhoidal Circular Stapler

| Code      | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|-----------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSSWBE32 | 32                  | 22                  | 32                | 3.8                     | 1.5                       | 0.28             | 1            |
| FCSSWBE33 | 33                  | 24.2                | 32                | 3.8                     | 1.5                       | 0.28             | 1            |
| FCSSWBE34 | 34                  | 24.2                | 34                | 3.8                     | 1.5                       | 0.28             | 1            |

## Haemorrhoidal Circular Stapler with Transparent Cartridge

| Code      | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|-----------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSSWAE32 | 32                  | 22                  | 32                | 3.8                     | 1.5                       | 0.28             | 1            |
| FCSSWAE33 | 33                  | 24.2                | 32                | 3.8                     | 1.5                       | 0.28             | 1            |
| FCSSWAE34 | 34                  | 24.2                | 34                | 3.8                     | 1.5                       | 0.28             | 1            |

## Haemorrhoidal Circular Stapler with Transparent Cartridge and Sound Feedback

| Code     | Outer Diameter (mm) | Inner Diameter (mm) | Number of Staples | Staple Height Open (mm) | Staple Height Closed (mm) | Staple Wire (mm) | Box Quantity |
|----------|---------------------|---------------------|-------------------|-------------------------|---------------------------|------------------|--------------|
| FCSSME32 | 32                  | 23.9                | 32                | 3.8                     | 1.5                       | 0.28             | 1            |
| FCSSME34 | 34                  | 25.9                | 34                | 3.8                     | 1.5                       | 0.28             | 1            |

# Hernia Fixation Device



| Code     | Shaft Length (mm) | Diameter (mm) | Staple Height (mm) | Tack Wire Diameter (mm) | Number of Tacks | Box Quantity |
|----------|-------------------|---------------|--------------------|-------------------------|-----------------|--------------|
| HRDN0530 | 355               | 5             | 3.8                | 4                       | 30              | 1            |



## Research & Development

Research driven and forward-thinking to produce the latest surgical stapling innovations.



## Sustainability

Committed to ensuring a sustainable healthcare ecosystem, Panther Healthcare is ISO 14001 certified.



## Panther Surgical Academy

Global programme of training courses with opportunities to learn from the most experienced of surgeons.



Panther Healthcare UK is operated by LaproSurge, the official UK distributor of Panther Healthcare's surgical stapling devices.

LaproSurge are a vastly experienced British medical device manufacturer and distributor founded in 1999 with an aim to provide better value laparoscopic devices to hospitals. Based in Watford, LaproSurge has grown into a company with a strong UK market presence and a large network of international distributors.

In the ever growing minimally invasive market, LaproSurge has consistently offered an extensive portfolio of high-quality single-use devices for laparoscopic surgery. LaproSurge and Panther Healthcare are strategically aligned to ensure that UK hospitals have available the latest generations of surgical stapling technology.

For all orders and enquiries  
contact LaproSurge at  
[info@laprosurge.com](mailto:info@laprosurge.com)



**UK Manufactured, Globally Distributed, Minimally Invasive.**

+44 (0) 1923 220 711

[info@laprosurge.com](mailto:info@laprosurge.com)

LaproSurge Ltd, Unit 5b, Fishers Industrial Estate, Wiggenhall Road,  
Watford, Hertfordshire, WD18 0FN, England

[laprosurge.com](https://laprosurge.com)

© 1999–2026 LaproSurge Ltd